

MICROSCOPE DAILY CLEANING CHECKLIST CANTON CITY PUBLIC HEALTH

MONTH: _____

YEAR: _____

After each day of use, indicate cleaning performed by placing a check mark in the appropriate column and initialing at the bottom of the daily column. Only check items that were done, not all items will need to be done every day. NOTE: When cleaning be very conservative in the application of liquids to all surfaces since moisture could damage components.

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
Turn off light																																
Use lens cleaning fluid on lens tissue to wipe lenses of dust debris or oil.																																
If body fluids are suspected of contacting scope, use bactericidal disinfectant to decontaminate.																																
Replace Dust Cover																																
INITIALS																																

Comments: