

**Canton City Health District  
Syphilis- Rapid Plasma Reagin (RPR) Quality Control / Quality Assurance Log**

**for year:** \_\_\_\_\_ **(to be performed each day of testing)**

Date	Temp. (73-84° F, 23-29° C)	Centriuge Secure (Y/N)	RPM (100 +/- 2)	Antigen: <b>B-D Macro-Vue</b> , lot #, expiration	Controls: <b>B-D Macro-Vue</b> , lot #, expiration	R / MR / Non	Pt. # / Titer	Titer Control / Titer	Anal.

RPM=Revolutions Per Minute, R/MR/Non= Reactive / Minimally-Reactive / Non-Reactive, Pt=Patient, Anal. = Analyst interpreting the results, API=American Proficiency Institute